
# seca 645

|           |  |            |
|-----------|--|------------|
| <b>Ⓚ</b>  | <b>Bedienungsanleitung und Garantieerklärung . . . . .</b>       | <b>3</b>   |
| <b>ⒼⒷ</b> | <b>Instruction manual and guarantee . . . . .</b>                | <b>61</b>  |
| <b>ⒻⓇ</b> | <b>Mode d'emploi et garantie . . . . .</b>                       | <b>117</b> |
| <b>ⒾⓉ</b> | <b>Manuale di istruzioni e garanzia . . . . .</b>                | <b>175</b> |
| <b>ⒺⓈ</b> | <b>Manual de instrucciones y garantia . . . . .</b>              | <b>233</b> |
| <b>Ⓟ</b>  | <b>Instruções de utilização e declaração de garantia . . . .</b> | <b>289</b> |
| <b>Ⓟ</b>  | <b>Instrukcja obsługi i gwarancja . . . . .</b>                  | <b>345</b> |
|           | <b>Declaration of conformity for USA and Canada . . . . .</b>    | <b>404</b> |

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# TABLE OF CONTENTS

|   |           |  |
|---|-----------|--|
| <b>1. Device description</b> .....                            | <b>63</b> |  |
| 1.1 Intended use .....  | 63        |  |
| 1.2 Description of function .....                             | 63        |  |
| 1.3 User qualification .....                                  | 63        |  |
| Assembly .....  | 63        |  |
| Operation .....   | 64        |  |
| <b>2. Safety information</b> .....                            | <b>64</b> |  |
| 2.1 Safety information in these<br>instructions for use ..... | 64        |  |
| 2.2 Basic safety instructions .....                           | 64        |  |
| Handling the instrument .....                                 | 64        |  |
| Preventing electric shock .....                               | 66        |  |
| Avoiding infections .....                                     | 66        |  |
| Avoiding injuries .....                                       | 67        |  |
| Preventing device damage .....                                | 67        |  |
| Handling measuring results .....                              | 68        |  |
| Handling packaging material .....                             | 69        |  |
| 2.3 Handling batteries and<br>battery packs .....             | 69        |  |
| <b>3. Overview</b> .....                                      | <b>71</b> |  |
| 3.1 Controls .....  | 71        |  |
| 3.2 Symbols in the display .....                              | 72        |  |
| 3.3 Marking on the device and<br>on the rating plate .....    | 73        |  |
| 3.4 Marking on the packaging .....                            | 74        |  |
| 3.5 Menu structure .....                                      | 75        |  |
| <b>4. Before you really get<br/>started...</b> .....          | <b>76</b> |  |
| 4.1 Scope of delivery .....                                   | 76        |  |
| Components .....  | 76        |  |
| Connecting elements .....                                     | 77        |  |
| 4.2 Assembling the device .....                               | 78        |  |
| Preparing the weighing<br>platform .....                      | 78        |  |
| Fitting the right-hand handrail<br>support .....              | 79        |  |
| Fitting the rear handrail<br>support .....                    | 80        |  |
| Fitting the handrail .....                                    | 81        |  |
| Fitting the display bracket .....                             | 82        |  |
| Fitting the display housing .....                             | 83        |  |
|   |           | Fitting the triangular cover<br>panels .....                         |
|   |           | 83   |
|   |           | Fitting the caster bar and<br>side bar .....                         |
|   |           | 84   |
|   |           | Connecting the display cable .....                                   |
|   |           | 85   |
|   |           | 4.3 Establishing power supply .....                                  |
|   |           | 86   |
|   |           | Connecting the battery block .....                                   |
|   |           | 86   |
|   |           | Connecting the power pack<br>and charging the battery<br>block ..... |
|   |           | 87   |
| <b>5. Operation</b> .....                                     | <b>88</b> |  |
| 5.1 Setting up the scale .....                                | 88        |  |
| Moving the scale .....  | 88        |  |
| Aligning the scale .....                                      | 88        |  |
| Turning the display head .....                                | 89        |  |
| 5.2 Weighing .....  | 90        |  |
| Switching on the scale .....                                  | 90        |  |
| Weighing a patient .....                                      | 90        |  |
| Taring off extra weight<br>(TARE) .....                       | 90        |  |
| Keeping the measurement<br>result in the display (HOLD) ..... | 91        |  |
| Determining and rating body<br>mass index (BMI) .....         | 92        |  |
| Sending measurement results<br>to wireless receivers .....    | 93        |  |
| Calculating and printing out<br>BMI automatically .....       | 93        |  |
| Automatic switching of<br>weighing range .....                | 94        |  |
| Switching off the scale .....                                 | 94        |  |
| 5.3 Other functions (menu) .....                              | 95        |  |
| Navigating in the menu .....                                  | 95        |  |
| Deleting saved values<br>automatically (ACLR) .....           | 96        |  |
| Permanently saving an<br>additional weight (Pt) .....         | 96        |  |
| Activating the auto hold<br>(AHOLD) function .....            | 98        |  |
| Activating beeps (bEEP) .....                                 | 98        |  |
| Setting filtering (FIL) .....                                 | 99        |  |
| Restoring factory settings<br>(rESET) .....                   | 99        |  |

|  |            |   |            |
|--|------------|---|------------|
| <b>6. The seca 360° wireless network</b> .....           | <b>100</b> | <b>8. Function check</b> .....                              | <b>108</b> |
| 6.1 Introduction .....                                   | 100        | <b>9. Troubleshooting</b> .....                             | <b>108</b> |
| seca wireless groups .....                               | 100        | <b>10. Servicing/Verification</b> .....                     | <b>110</b> |
| Channels .....   | 101        | 10.1 Information on servicing and verification .....        | 110        |
| Device detection .....                                   | 101        | 10.2 Checking the content of the verification counter ..... | 111        |
| 6.2 Operating the scale in a wireless group (menu) ..... | 102        | <b>11. Technical data</b> .....                             | <b>112</b> |
| Activating the wireless module (SYS) .....               | 102        | 11.1 General technical data .....                           | 112        |
| Setting up a wireless group (Lrn) .....                  | 102        | 11.2 Technical data, weighing .....                         | 113        |
| Activating automatic send (ASEnd) .....                  | 104        | <b>12. Optional accessories</b> .....                       | <b>113</b> |
| Selecting print option (APrt) .....                      | 105        | <b>13. Spare parts</b> .....                                | <b>114</b> |
| Setting the time (TIME) .....                            | 105        | <b>14. Disposal</b> .....                                   | <b>114</b> |
| <b>7. Hygiene treatment</b> .....                        | <b>106</b> | 14.1 Device .....   | 114        |
| 7.1 Cleaning .....                                       | 107        | 14.2 Batteries and rechargeable batteries .....             | 114        |
| 7.2 Disinfecting .....                                   | 107        | <b>15. Warranty</b> .....                                   | <b>115</b> |
| 7.3 Sterilizing .....                                    | 107        | <b>16. Declaration of conformity</b> ...                    | <b>115</b> |

# 1. DEVICE DESCRIPTION

## 1.1 Intended use

---

The **seca 645** electronic scale with standing aid is mainly used in hospitals, doctors' surgeries and in-patient care facilities in accordance with national regulations.

The scale is for conventional determination of weight and establishment of general state of nutrition; it assists the physician supervising treatment in making a diagnosis or deciding on a course of treatment.

To make an accurate diagnosis, however, other specific examinations have to be ordered by the physician and their results taken into account, in addition to determining a weight value.

## 1.2 Description of function

---

In addition to conventional determination of weight, the **seca 645** provides a function for determining body mass index. To this end, height is entered on the keypad and the body mass index associated with the weight value is calculated automatically. Length measuring devices from the **seca 360° wireless** system can transmit height to the **seca 645** wirelessly.

Via the **seca 360° wireless** network, measuring results can be transmitted wirelessly to a seca wireless printer or to a PC equipped with **seca analytics** PC software and the seca USB wireless adapter.

The **seca 645** can be moved on casters.

Use the scale only for the purpose named in the section entitled "Intended use" on page 63.

## 1.3 User qualification

---

### Assembly

Devices that are shipped partially assembled may only be mounted by sufficiently qualified persons such as specialist dealers, hospital technicians or seca service technicians.

**Operation** The device may only be operated by healthcare professionals.

## 2. SAFETY INFORMATION

### 2.1 Safety information in these instructions for use

---

**DANGER!**

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries will occur.

**WARNING!**

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries may result.

**CAUTION!**

Used to identify a hazardous situation. If you fail to take note of this information, minor to moderate injuries may result.

**NOTICE!**

Used to identify possible incorrect usage of the device. If you fail to take note of this information, you may damage the device, or the measured results may be incorrect.

**NOTE**

Includes additional information about use of the device.

### 2.2 Basic safety instructions

---

**Handling the instrument**

- ▶ Please take note of the information in these instructions for use.
- ▶ Keep the instructions for use in a safe place. The instructions for use are a component of the device and must be available at all times.

**DANGER!****Risk of explosion**

Do not use the device in an environment in which one of the following gases has accumulated:

- oxygen
- flammable anesthetics
- other flammable substances/air mixtures

**CAUTION!****Patient hazard, damage to device**

- ▶ Additional devices which are connected to electrical medical devices must provide evidence of compliance with the relevant IEC or ISO standards (e.g. IEC 60950 for data-processing devices). Furthermore, all configurations must comply with the requirements of standards for medical systems (see IEC 60601-1-1 or Section 16 of the 3rd edition of IEC 60601-1 respectively). Anyone connecting additional devices to electrical medical devices is considered a system configurator and is therefore responsible for ensuring that the system complies with the requirements of standards for systems. This also applies to additional devices recommended by seca. Your attention is drawn to the fact that local laws take precedence over the above-mentioned requirements of standards. In the event of any queries, please contact your local specialist dealer or Technical Service.
- ▶ Have servicing and verification carried out regularly as described in the relevant section of this document.
- ▶ Technical modifications may not be made to the device. The device does not contain any parts for servicing by the user. Only have servicing and repairs performed by an authorized seca Service partner. You can find service partners in your area at [www.seca.com](http://www.seca.com) or by sending an e-mail to [service@seca.com](mailto:service@seca.com).
- ▶ Only use original seca accessories and spare parts, otherwise seca will not grant any warranty.

### **CAUTION!**

#### **Patient hazard, malfunction**

- ▶ Keep other electrical medical devices, e.g. high-frequency surgical devices, a minimum distance of approx. 1 meter away to prevent incorrect measurements or wireless transmission interference.
- ▶ Keep HF devices such as cell phones a minimum distance of approx. 1 meter away to prevent incorrect measurements or wireless transmission interference.
- ▶ The actual transmission output of HF equipment may require minimum distances of more than 1 meter. Details can be found at [www.seca.com](http://www.seca.com).

## **Preventing electric shock**

### **WARNING!**

#### **Electric shock**

- ▶ Set up the device so that the power supply socket is easy to reach and the device can be disconnected from the power supply quickly.
- ▶ Ensure that your local power supply matches the information on the power supply unit.
- ▶ Do not touch the power supply unit with wet hands.
- ▶ Do not use extension cables or power strips.
- ▶ Make sure that cables are not pinched or damaged by sharp edges.
- ▶ Make sure that cables do not come into contact with hot objects.
- ▶ Do not operate the device at an altitude of more than 3000 m above sea level.

## **Avoiding infections**

### **WARNING!**

#### **Risk of infection**

- ▶ Before and after every measurement, wash your hands to reduce the risk of cross-contamination and nosocomial infections.
- ▶ Hygienically reprocess the scale regularly as described in the respective section in this document.
- ▶ Make sure that the patient has no infectious diseases.

- ▶ Make sure that the patient has no open wounds or infectious skin alterations, which may come into contact with the device.

## Avoiding injuries

### **WARNING!** **Injury from falls**

- ▶ Ensure that the device is positioned firmly and level.
- ▶ Route connecting cables (if present) in such a way that neither user nor patient can trip over them.
- ▶ Make sure that the patient does not step onto and off the weighing platform right at the edges.
- ▶ Make sure that the patient steps onto and off the weighing platform slowly and safely.

### **WARNING!** **Risk of slipping**

- ▶ Ensure that the weighing platform is dry before the patient steps onto it.
- ▶ Ensure that the patients feet are dry before he or she steps onto the weighing platform.
- ▶ Make sure that the patient steps onto and off the weighing platform slowly and safely.

## Preventing device damage

### **NOTICE!** **Damage to device**

- ▶ Ensure that no liquids enter the device. They can damage the electronics.
- ▶ Switch off the device before disconnecting the power supply unit from the mains socket.
- ▶ For devices with power supply operation: Disconnect the power supply unit from the mains socket if you intend to not use the device for a longer period of time. Only this way it can be ensured that the device is de-energized.
- ▶ For devices with battery or rechargeable battery operation: If you do not use the device for an extended period of time, remove batteries or rechargeable batteries. Only this way it can be ensured that the device is de-energized.
- ▶ Make sure not to drop the device.

- ▶ Do not expose the device to any impacts or vibrations.
- ▶ Perform function controls regularly as described in the relevant section in this document. Do not operate the device if it is damaged or not working properly.
- ▶ Ensure that there is no heat source in the immediate vicinity. Do not expose to direct sunlight. The excessive temperature could damage the electronics.
- ▶ Avoid rapid temperature fluctuations. When the device is transported so that a temperature difference of more than 20 °C occurs, it must stay turned off for at least 2 hours before it can be turned on again. Otherwise, condensation water will form which can damage the electronics.
- ▶ Use the device only in the ambient conditions outlined in "Intended use".
- ▶ Store the device only in the storage conditions outlined in "Intended use".
- ▶ Use only chlorine and alcohol-free disinfectants which are explicitly suitable for acrylic sheet and other sensitive surfaces (active ingredient: quaternary ammonium compounds, for example).
- ▶ Do not use aggressive or abrasive cleaning agents.
- ▶ Do not use organic solvents (e.g. white spirit or petroleum spirit).

## Handling measuring results

### **WARNING!** **Patient hazard**

This device is **no** diagnostic device. It simply assists the treating physician in establishing a diagnosis.

- ▶ In order to make a precise diagnosis and initiate therapeutic measures, besides determination of the weight, further targeted examinations must be set up by the physician, and their results must be considered.
- ▶ The responsibility for diagnosis and treatment lies with the treating physician.

**NOTICE!****Inconsistent measuring results**

- ▶ Before you electronically save measurement values determined using this device and use them further (e.g. in seca PC software or in a hospital information system), make sure that the measurement values are plausible.
- ▶ If measurement values are transmitted to seca PC software or a hospital information system, make sure prior to further use that the measurement values are plausible and are assigned to the correct patient.

**Handling packaging material****WARNING!****Risk of suffocation**

Packaging material made of plastic foil (bags) is a choking hazard.

- ▶ Keep packaging material out of reach of children.
- ▶ In the event that the original packing material may not be available anymore, only use plastic bags with security holes in order to reduce the risk of suffocation. Use recyclable materials if possible.

**NOTE**

Keep the original packing material for future use (e.g. returning for servicing).

## 2.3 Handling batteries and battery packs

---

**WARNING!****Personal injury as a result of improper handling**

Batteries and rechargeable batteries contain harmful substances which may explode if not handled properly.

- ▶ Do not try to recharge batteries.
- ▶ Do not expose (rechargeable) batteries to heat.
- ▶ Do not burn (rechargeable) batteries.
- ▶ If acid is leaking out, avoid contact with the skin, eyes and mucous membranes. Rinse affected areas with plenty of clean water and seek medical help at once.


## **NOTICE!**

### **Damage to device and malfunctions with improper handling**

- ▶ Only use the type of (rechargeable) battery specified in this document.
- ▶ When replacing (rechargeable) batteries, always replace a complete set at a time.
- ▶ Do not short-circuit (rechargeable) batteries.
- ▶ If you do not use the device for a long period of time, remove the batteries (incl. rechargeable batteries). This prevents acid from leaking into the device.
- ▶ If acid leaked into the device, discontinue use. Have the device checked by an authorised service partner and repaired if necessary.

## 3. OVERVIEW

### 3.1 Controls

| No. | Control   | Function  |
|-----|---|---|
| 1   | Display housing   | Central control and display element   |
| 2   |  | Switch scale on and off   |
| 3   |  | Arrow key <ul style="list-style-type: none"> <li>• During weighing:               <ul style="list-style-type: none"> <li>- press briefly: activate hold function</li> <li>- press and hold: activate tare function</li> </ul> </li> <li>• In the menu:               <ul style="list-style-type: none"> <li>- select submenu, select menu item</li> <li>- increase value</li> </ul> </li> </ul> |

| No. | Control   | Function   |
|-----|---|--|
| 4   |  | <p>Arrow key</p> <ul style="list-style-type: none"> <li>• During weighing: <ul style="list-style-type: none"> <li>- press briefly: activate BMI function</li> <li>- press and hold: call up menu</li> </ul> </li> <li>• In the menu: <ul style="list-style-type: none"> <li>- select submenu, select menu item</li> <li>- reduce value</li> </ul> </li> </ul>  |
| 5   |  | <p>Enter key</p> <ul style="list-style-type: none"> <li>• During weighing (if wireless network is set up): <ul style="list-style-type: none"> <li>- press briefly: send measuring result to devices ready to receive it (PC with USB wireless adapter)</li> <li>- press and hold: print out measuring result (wireless printer)</li> </ul> </li> <li>• In the menu: <ul style="list-style-type: none"> <li>- confirm selected menu item</li> <li>- save set value</li> </ul> </li> </ul> |
| 6   | Display   | Display element for measuring results and to configure the device  |
| 7   | Transport casters   | The scale can be moved on these casters.   |
| 8   | Power supply connection   | Serves to connect the power pack supplied  |
| 9   | Spirit level  | Shows whether the device is horizontal   |
| 10  | Battery compartment   | Contains the battery block   |
| 11  | Foot screw  | 4 pcs, for precise alignment   |

## 3.2 Symbols in the display


|   | Symbol  | Meaning   |
|---|---|---|
| A |  | Operation with power pack   |
| B |  | Non-calibratable function active                                    |
| C |  | Memory location currently in use                                    |
| D |  | Weighing range currently in use<br>see "Technical data" on page 112 |

### 3.3 Marking on the device and on the rating plate

| Text/symbol   | Meaning  |
|---|--|
|    | Name and address of manufacturer, date of manufacture  |
| <b>REF</b>  | Model number   |
| <b>SN</b>   | Serial number, consecutive   |
| <b>GAL</b>  | Value in $m/s^2$ (verified models) <ul style="list-style-type: none"> <li>• Gravitational acceleration on earth</li> <li>• Depends on the intended location</li> </ul>   |
| <b>ProdID</b>   | Product identification number, consecutive   |
| <b>Approval Type</b>  | Type designation of design approval  |
|    | Follow instructions for use  |
|    | Electrical medical device, type B  |
| e   | Value in units of measurement used to classify and verify a scale  |
| d   | Value in units of measurement which gives the difference between two consecutive display values  |
| $\rightarrow x \leftarrow$  | Active weighing range  |
|    | Class III scale to directive 2014/31/EU  |
|  | Device complies with EU standards and directives. <ul style="list-style-type: none"> <li>• <b>M</b>: Conformity label according to Directive 2014/31/EU governing non-automatic weighing instruments (verified models)</li> <li>• <b>18</b>: (Example: 2018) Year in which conformity verification was performed and the CE label was applied (verified model)</li> <li>• <b>0102</b>: Notified body metrology (verified models)</li> <li>• <b>0123</b>: Notified body medical products</li> </ul> |
| <b>FC</b>   | Symbol of the US Federal Communications Commission (FCC)   |
| FCC ID  | Device license number from the Federal Communications Commission (FCC)   |
| IC  | Device license number from Industry Canada   |


| Text/symbol   | Meaning   |
|---|---|
|  <p>x-y V <math>\equiv</math><br/>max. xxx mA<br/>use compatible<br/>seca adapter only</p> | <p>Rating plate on the power supply connection socket</p> <ul style="list-style-type: none"> <li>• <b>x-y V</b>: required supply voltage</li> <li>• <b>max xxx mA</b>: maximum current consumption</li> <li>• <math>\text{--}\ominus\text{+}</math> : note polarity of device plug</li> <li>• <math>\equiv</math> : operate device with direct current</li> </ul> |
|    | <p>Do not dispose of device with household waste</p>  |

### 3.4 Marking on the packaging

|   |  |
|---|--|
|    | <p>Protect from moisture</p>   |
|    | <p>Arrows indicate top of product<br/>Transport and store in an upright position</p> |
|    | <p>Fragile<br/>Do not throw or drop</p>  |
|    | <p>Permitted min. and max. temperature for transport and storage</p>                 |
|    | <p>Permitted min. and max. moisture for transport and storage</p>                    |
|    | <p>Not sterile</p>   |
|  | <p>Do not reuse</p>  |
|  | <p>Open packaging here</p>   |
|  | <p>Packaging material can be disposed of through recycling programs</p>              |

### 3.5 Menu structure

Other functions are available to you in the menu for the device. This enables you to configure the device perfectly to suit your needs (details from page 95).

**\*seca 360° wireless network:**

**Group (ID)**  
 Max. three seca wireless groups: 0, 1, 2  
**Maximum configuration per wireless group**

- 1 baby scale
- 1 personal scale
- 1 length measuring rod
- 1 wireless printer
- 1 PC with USB wireless adapter

**Channel (C1, C2, C3)**

- Three channels per wireless group (total 9 channels)
- Channel numbers: 0 - 99
- Only use each channel number once
- Recommended spacing: 30

**Example configuration**

- Group 0: C1 0, C2 30, C3 60
- Group 1: C1 10, C2 40, C3 70
- Group 2: C1 20, C2 50, C3 80

(Note: no spaces in display)


**Registered devices (MO)**

- 1: Personal scale
- 2: Length measuring rod
- 3: Wireless printer
- 4: PC with USB wireless adapter
- 7: Baby scale

**Call up menu:**


**Navigation:**





## 4. BEFORE YOU REALLY GET STARTED...

### 4.1 Scope of delivery

| No. | Component  | Pcs. |
|-----|--|------|
| A   | Rear handrail support with display head                                  | 1    |
| B   | Right-hand railing support with opening for power pack connection        | 1    |
| C   | Weighing platform  | 1    |
| D   | Handrail   | 1    |
| E   | Caster bar   | 1    |
| F   | Side bar with opening for power pack connection                          | 1    |
| G   | Display bracket  | 1    |
| H   | Power pack with adapters (depending on model: Euro connector power pack) | 1    |
|     | Instructions for use, not shown  | 1    |

## Connecting elements


| No.      | Component                              | Pcs. |
|----------|--|------|
| <b>a</b> | Hex socket head screw M 5 x 12         | 6    |
| <b>b</b> | Pan-head tapping screw 6.3 x 38        | 1    |
| <b>c</b> | Raised countersunk head screw M 6 x 40 | 2    |
| <b>d</b> | Hex socket head screw M 6 x 50         | 1    |
| <b>e</b> | Pan-head tapping screw 2.9 x 9.5       | 3    |
| <b>f</b> | Pan-head tapping screw 3.5 x 9.5       | 16   |
|          | Allen key, 4 mm, not shown             | 1    |
|          | Phillips screwdriver size 1, not shown | 1    |
|          | Phillips screwdriver size 3, not shown | 1    |

## 4.2 Assembling the device

We recommend that two people assemble the scale, as large individual parts have to be aligned and then screwed together.

### Preparing the weighing platform

1. Place the weighing platform on a flat, stable surface.
2. Align the weighing platform as shown in the diagram above.


The weighing platform is in the correct position for fitting the handrail supports as shown in the following sections.

## Fitting the right-hand handrail support


You need the following connecting elements for this assembly step:

3 x hex socket head screws, M 5 x 12

1. Hold the right-hand handrail support in position on the weighing platform as shown in the diagram above.
2. Attach the handrail support to the weighing platform with two hex socket head screws.


3. Attach the triangular brace to the weighing platform with one hex socket head screw.

## Fitting the rear handrail support


You need the following connecting elements for this assembly step:


3 x hex socket head screws, M 5 x 12

1 x pan-head tapping screw 6.3 x 38

1. Hold the rear handrail support in position on the weighing platform as shown in the diagram above.
2. Attach the rear handrail support to the weighing platform with two hex socket head screws.


3. Attach the triangular brace to the weighing platform with one hex socket head screw.
4. Join the two handrail supports using a pan-head tapping screw as shown in the diagram above.


**Fitting the handrail**

You need the following connecting elements for this assembly step:


2 x raised countersunk head screws, M 6 x 40


1. Align the handrail so that the large bores point downwards.
2. Put the handrail on the three columns.
3. Fix the handrail in position on the two outer columns using one raised countersunk head screw for each.

## Fitting the display bracket

You need the following connecting element for this assembly step:

1 x cylinder bolt, M 6 x 50


1. Put the bracket holder in the bore above the central column.
2. Turn the display bracket until it engages audibly.
3. Fix the display bracket in position using the cylinder bolt.

## Fitting the display housing

You need the following connecting elements for this assembly step:

3 x pan-head tapping screws 2.9 x 9.5


1. Remove the transport securing device of the display housing.
2. Position the display housing on the display bracket in such a way that the position markings line up with one another.
3. Fix the display housing to the display bracket with the pan-head tapping screws.

## Fitting the triangular cover panels

You need the following connecting elements for this assembly step:

6 x pan-head tapping screws 3.5 x 9.5


1. Place one panel on the triangular brace for the rear handrail support.
2. Fix the panel to the brace with three pan-head tapping screws.
3. Repeat steps 1. and 2. for the brace for the right-hand handrail support.

## Fitting the caster bar and side bar


You need the following connecting elements for this assembly step:

10 x pan-head tapping screws 3.5 x 9.5

1. Place the caster bar on the rear handrail support as shown in the diagram above.
2. Fix the caster bar to the rear handrail support with 5 pan-head tapping screws.

3. Repeat steps 1. and 2. to fit the side bar to the right-hand handrail support.

## Connecting the display cable

### ATTENTION!

#### Malfunction due to installation errors


If cables are strained during fitting, e. g. with sharp bends or kinked connectors, this may result in faulty displays and failure of the display.

- ▶ Route all cables to prevent sharp bends and kinked connectors.
- ▶ Provide strain relief by routing all cables in the relevant holders.

1. Lift the scale so that the underside of the weighing platform is accessible.

### NOTE

You can rest the scale on the handrail for this assembly work.


2. Pass the display cable along under the weighing platform to the electronics box as shown in the diagram.
3. Connect the display cable to the electronics box.
4. Push the cable clip attached to the display cable at the factory into the bore on the scale frame until you feel it engage.
5. Press the display cable into the cable clips on the scale frame until you feel them engage.
6. Lower the scale.

## 4.3 Establishing power supply

The scale is powered by a battery block or a power pack (both included in scope of supply). Depending on model, the scope of supply will include either a power pack with plug adapters or a power pack with a fixed Euro connector.

### Connecting the battery block

1. Undo the screws of the battery compartment.
2. Remove the lid of the battery compartment.
3. Remove the battery block from the battery compartment.
4. Remove the connecting cable from the battery compartment.
5. Connect the connecting cable to the battery block.
6. Insert the battery block in the battery compartment.
7. Screw the battery compartment lid onto the battery compartment.

## Connecting the power pack and charging the battery block



### WARNING!

Using the wrong power packs may cause bodily injury or damage to the device

Conventional power packs may deliver a higher voltage than is indicated on them. The scale may overheat, catch fire, melt or short-circuit.

- ▶ Use only original seca plug-in power packs with a 9 V or controlled 12 V output voltage.

1. If required, insert the power supply connector necessary for your power supply into the power pack.
2. Insert the connector of the power pack into the connecting socket of the scale.
3. Plug the power pack into a power supply socket.
4. Leave the scale connected to the power supply for at least 24 hours when charging for the first time in order to charge the battery block fully.

## 5. OPERATION


### **CAUTION!** **Bodily injury**

Perform a function check as described in the section entitled "Function check" on page 108 before using the device each time.

### 5.1 Setting up the scale

---

#### Moving the scale

1. Tilt the device until it can move freely on the casters.
2. Keep the device in this position and move it to the desired set-up or storage location.

#### Aligning the scale



### **ATTENTION!**

#### **Incorrect measurement due to force shunts**

If the scale and housing is in contact with something, e.g. a towel, weight will not be measured correctly.

- ▶ Set the scale up so that only its foot screws are in contact with the floor.
1. Place the scale on firm, level surface.

### Turning the display head


2. Undo the knurled wheels.
3. Level the device by turning the foot screws.

The air bubble in the spirit level must be located in the exact center of the circle.

4. Tighten the knurled wheels in the direction of the arrow.

The foot screws are secured against being adjusted.

The display housing of the device swivels. This allows you to adjust it to suit every application situation perfectly.

- ▶ Swivel the display housing so that you can operate and read it comfortably.

## 5.2 Weighing

### CAUTION!

#### Patients can be injured if they fall

People with restricted mobility may fall if they try to step onto the scale.

- ▶ Support people with restricted mobility when they step onto the scale.

### Switching on the scale

- ▶ Press the Start key.  
All the elements of the display are shown briefly, then **SECA** appears in the display.  
The scale is operational when **0.0** appears in the display.

### Weighing a patient

1. Ensure that the scale has no load.
2. Ask the patient to step onto the scale.
3. Ask the patient to keep still.
4. Read off the measurement result.

### Taring off extra weight (TARE)

Using the TARE function, you can prevent any extra weight (e.g. a towel or a support on the weighing platform) from affecting the weighing result.

### ATTENTION!

#### Incorrect measurement due to force shunts

If an extra weight, e.g. a large towel, contacts the surface on which the scale is placed, weight will not be measured correctly.

- ▶ Make sure that extra weights are only placed on the scale's weighing surface.

1. Switch on the scale.
2. Place the extra weight on the scale.
3. Hold the **(hold/tare)** arrow key down until the message "NET" appears in the display.
4. Wait until the display stops flashing and **0.0** appears instead.
5. Weigh the patient as described in the section entitled "Weighing a patient".

6. Read off the measurement result.  
The additional weight is automatically deducted.
7. To deactivate the TARE function, press the **(hold/tare)** arrow key until the message "NET" is no longer displayed or turn off the scale.

#### NOTE

The maximum weight which can be displayed is reduced by the weight of the objects already placed on the scale.

### Keeping the measurement result in the display (HOLD)

When you activate the HOLD function, the weight value continues to be displayed after the weight has been removed from the scale. As a result, you can attend to the patient before recording the weight.

1. Ensure that the scale has no load.
2. Switch on the scale.
3. Weigh the patient as described in the section entitled "Weighing a patient".
4. Press the **(hold/tare)** arrow key briefly.

The display flashes until a stable weight is measured. The weight is then continuously displayed. The  $\triangle$  symbol (non-calibratable function) and the message "HOLD" are displayed.


5. To deactivate the HOLD function, briefly press the **(hold/tare)** arrow key.  
The  $\triangle$  symbol and the "HOLD" message are no longer displayed.

#### NOTE

If the auto hold function is activated, the weight value is automatically displayed permanently until the scale switches itself off or is switched off, see "Activating the auto hold (AHOLD) function" on page 98.

## Determining and rating body mass index (BMI)

Body mass index relates height and weight to one another and consequently allows more accurate information than ideal Broca weight, for example. A tolerance range is quoted which is considered ideal for health.

The device has three memories for height. You can enter and save the height of specific patients. Alternatively, you can save different start values and thus set a patient's actual height more quickly.

1. Ensure that the scale has no load.
2. Switch on the scale.
3. Press the **(bmi/menu)** arrow key briefly. The message "BMI" appears. The last memory used is displayed (in this case: memory 2).
4. You can adopt the memory displayed or select a different memory using the arrow keys.

5. Confirm your setting with Enter **(send/print)**. Arrows flash in the display.

The last height saved to the selected memory is displayed.

6. You can accept the height displayed or use the arrow keys to set a different height.

7. Confirm your setting with Enter **(send/print)**. The height entered is saved and will be available again for the next BMI calculation.

### NOTE

Make a note of the memory number so that you can call up this height for another BMI calculation.

8. Weigh the patient as described in the section entitled "Weighing a patient". The patient's BMI is automatically calculated and displayed.
9. Read off the BMI and compare it with the categories quoted further down.

10. To deactivate the BMI function, briefly press the Enter key (**send/print**).

| BMI  | Rating   |
|--|--|
| <b>Below 18.5</b>                            | Patient is underweight. There could be anorexic tendencies. Weight should be increased to improve well-being and performance. If in doubt, seek the advice of a consultant.  |
| <b>Between 18.5 and 24.9</b>                 | The patient's weight is normal.  |
| <b>Between 25 and 30 (incipient obesity)</b> | The patient is slightly to moderately overweight. He or she should lose weight if already suffering from a disease such as diabetes, high blood pressure, gout or disorders of fat metabolism.   |
| <b>Over 30</b>                               | Weight needs to be lost urgently. The metabolism, the circulatory system and the bones are all under strain. A strict diet, lots of physical activity and behavioral training are all recommended. If in doubt, seek the advice of a consultant. |

### Sending measurement results to wireless receivers

If the scale is integrated in a **seca 360° wireless** network, you can send the measurement results to devices ready to receive them (e.g. wireless printer, PC with USB wireless adapter) at the touch of a button.

- ▶ Press the Enter key (**send/print**).
  - Press key briefly: send measurement results to all devices ready to receive them
  - Press and hold key: print out measurement result on wireless printer

### Calculating and printing out BMI automatically

If you use this scale together with a wireless printer and a length measuring rod from the **seca 360° wireless** system, you can have BMI calculated and printed out automatically.

#### NOTE

The prerequisite for this function is that the devices are registered together in one wireless group (see “The seca 360° wireless network” on page 100).

1. Perform the height measurement.
2. Briefly press the Enter key (**send/print**) of the length measuring rod.  
The measured value is sent to the wireless printer but not printed.
3. Perform the weighing operation.

4. Press and hold the Enter key (**send/print**) of the scale.

The measured value is sent to the wireless printer.

BMI is calculated.

Height, weight and BMI are printed out.

### **Automatic switching of weighing range**

The scale has two weighing ranges. In weighing range 1 (→1←), capacity is reduced, but the increments in the weight display are finer. In weighing range 2 (→2←), you can use the maximum capacity of the scale.

After the scale is switched on, weighing range 1 is active. If a particular weight is exceeded, the scale automatically switches to weighing range 2.

To switch back to weighing range 1, proceed as follows.

- ▶ Completely remove the load from the scale.  
Weighing range 1 is active again.

### **Switching off the scale**


- ▶ Press the Start key.

#### **NOTE**

When in battery mode, the scale will switch off automatically after a short period if no load is placed on it.

## 5.3 Other functions (menu)

Other functions are available to you in the menu for the scale. This enables you to configure the scale to suit your conditions of use perfectly.


\* The description of menu item “rF” can be found in the section entitled “Operating the scale in a wireless group (menu)” on page 102.

### Navigating in the menu

1. Switch on the scale.
2. Keep arrow key (**bmi/menu**) pressed until the menu is called up.  
The last menu item selected (in this case: auto hold “AHOLD”) appears in the display.
3. Press one of the arrow keys until the desired menu item appears in the display (in this case: filtering “FIL”).
4. Confirm your selection with Enter (**send/print**).  
The current setting for the menu item or a submenu are displayed (in this case: level “0”).

- To amend the setting or call up a different sub-menu, press one of the arrow keys until the desired setting (in this case: level “2”) is displayed.
- Confirm the setting with Enter (**send/print**). You exit the menu automatically.
- To make more settings, call up the menu again and proceed as described.

#### NOTE

If no key is pressed for about 24 seconds, you exit the menu automatically.

- Press the Start key.  
The settings will be saved.  
The device switches off.

### Deleting saved values automatically (ACLR)

To prevent current measurement results remaining in the memory of the device and leading to an incorrect BMI calculation, you can set the scale so that measurement results are automatically deleted after 5 minutes.

#### NOTE


This function is activated at the factory on some models. You can deactivate the function if desired.

- Select the item “ACLR” from the menu.
- Confirm your selection.
- Select the desired setting:
  - On
  - Off
- Confirm your selection.  
You exit the menu automatically.

### Permanently saving an additional weight (Pt)

Using the Pre-tare function, you can permanently save an additional weight and automatically deduct it from a measurement result. For example, you can save an all-inclusive weight for shoes and clothing and then always deduct this from the measurement result if a patient is weighed fully-clothed.

The device has three memories for weights. You can save various weights and call them up individually depending on the initial situation, so that they are automatically deducted from the measurement result.

1. Select the item "Pt" from the menu.  
The message "Pt" appears.
2. Confirm your selection.  
The last memory used is displayed.
3. You can adopt the memory displayed or select a different memory.
4. Confirm your selection.  
Arrows flash in the display.  
The additional weight saved to the selected memory is displayed.
5. You can adopt the saved value or amend it using the arrow keys.

#### NOTE

When you enter the value "0", the function is switched off. The message "Pt" is no longer displayed.

6. Confirm your selection.
7. Ask the patient to position him or herself on the scale.  
The patient's weight is displayed.  
The saved additional weight has been deducted automatically.
8. To deactivate the function, select the item "Pt" from the menu again.
9. Confirm your selection.  
The function is deactivated.  
You exit the menu automatically.

#### NOTE

When you switch the scale off, the function is switched off. The message "Pt" is no longer displayed when you switch the scale on again.

## Activating the auto hold (AHOLd) function

When you activate the auto hold function, the measurement result for every weighing operation continues to be displayed after the weight has been removed from the scale. It is then no longer necessary to activate the Hold function manually for each individual weighing operation.

### NOTE

This function is activated at the factory on some models. You can deactivate the function if desired.

AHOLd

On

1. Select the item "AHOLd" from the menu.
2. Confirm your selection.  
The current setting is displayed.
3. Select the desired setting:
  - On
  - Off
4. Confirm your selection.  
You exit the menu automatically.

## Activating beeps (bEEP)

You can set whether a beep is heard every time a key is pressed and when a stable weight value is achieved. The latter is significant for the Hold/auto hold function.

### NOTE

The "Beep when weight is stable" function is activated at the factory. You can deactivate this function if desired.

bEEP

PrESS

On

1. Select the item "bEEP" from the menu.
2. Confirm your selection.
3. Select a menu item.
  - PrESS: beep when key is pressed
  - Hold: beep when weight value is stable
4. Confirm your selection.  
The current setting is displayed.
5. Select the desired setting:
  - On
  - Off
6. Confirm your selection.  
You exit the menu automatically.
7. If you also want to activate beeps for the second function, repeat the process.

**Setting filtering (FIL)**

Filtering (FIL) allows you to reduce interference (e.g. as a result of patient movement) when determining weight.

FIL

FIL 0

FIL 2

1. Select the item "FIL" from the menu.
2. Confirm your selection.

The current setting is displayed.

3. Select a filtering level.
  - 0: no filtering
  - 1: moderate filtering
  - 2: high filtering
4. Confirm your selection.  
You exit the menu automatically.

**Restoring factory settings (rESET)**

You can restore the factory settings for the functions below.

| Function                         | Factory setting |
|----------------------------------|-----------------|
| Auto hold (AHOLd)                | Model-dependent |
| Beep (PrESS)                     | Off             |
| Beep (Hold)                      | On              |
| Filtering (FIL)                  | 0               |
| Autoclear (Aclear)               | Model-dependent |
| Pre-tare (Pt)                    | 0 kg            |
| Height for Body Mass Index (BMI) | 170 cm          |
| Wireless module (SYS)            | Off             |
| Autosend (ASEnd)                 | Off             |
| Autoprint (APrt)                 | Off             |

**NOTE**

The wireless module is switched off when factory settings are restored. Information about existing wireless groups is retained. Wireless groups do not have to be set up again.

rESET

1. Select the item "rESET" from the menu.
2. Confirm your selection.  
You exit the menu automatically.
3. Switch off the scale.  
Factory settings are restored and are available when the scale is switched back on.

## 6. THE SECA 360° WIRELESS NETWORK

### 6.1 Introduction

This device is equipped with a wireless module. The wireless module allows wireless transmission of measurement results for evaluation and documentation. Data can be transmitted to the following devices:

- seca wireless printer
- PC with seca USB wireless adapter

#### **seca wireless groups**


The **seca 360° wireless** network operates with wireless groups. A wireless group is a virtual group of transmitters and receivers. If several transmitters and receivers of the same type are going to be operated, up to 3 wireless groups (0, 1, 2) can be set up with this device.

Setting up several wireless groups ensures the reliable transmission of correctly-addressed measured values if you wish to use several examination rooms each with comparable equipment.

The maximum distance between transmitters and receivers is approx. 10 meters. Certain local circumstances such as the thickness and type of walls, may reduce the range.


The following combination of devices is possible for each wireless group:

- 1 baby scale
- 1 personal scale
- 1 length measuring rod
- 1 seca wireless printer
- 1 PC with seca USB wireless adapter

## Channels

Within a wireless group, devices communicate with each other on three channels (C1, C2, C3). This guarantees that data are transmitted reliably and without interference.

If you set up a wireless group with this scale, the device suggests three channels which will ensure optimal data transmission. We recommend adopting the channel numbers suggested.

You can also manually select the channel numbers (0 to 99), for example if you wish to set up several wireless groups.

To ensure interference-free data transmission, the channels must be spaced sufficiently far apart. We recommend a spacing of min. 30 between channel numbers. Each channel number may only be used for one channel.

Example configuration for channel numbers when setting up 3 wireless groups within a medical practice:

- wireless group 0: C1=\_0, C2=30, C3=60
- wireless group 1: C1=10, C2=40, C3=70
- wireless group 2: C1=20, C2=50, C3=80


## Device detection

If you set up a wireless group with the scale, it searches for other active devices from the **seca 360° wireless** system. The devices detected are shown in the display of the scale in the form of modules (e.g. "MO 3"). The numbers have the following meaning:

- 1: personal scale
- 2: length measuring rod
- 3: wireless printer
- 4: PC with seca USB wireless adapter
- 7: baby scale
- 5, 6 and 8-12: reserved for system expansion

## 6.2 Operating the scale in a wireless group (menu)

All the functions you need to operate the device in a seca wireless group can be found in the “rF” submenu. Information for navigating in the menu can be found on page 95.

### Activating the wireless module (SYS)

The device is supplied with the wireless module deactivated. You have to activate it before you can set up a wireless group.


1. Switch on the device.
2. Select the “SYS” menu item from the “rF” submenu.
3. Confirm your selection.
4. Select the “On” setting.
5. Confirm your selection.  
You exit the menu automatically.
6. If you do not wish to make any more settings, press the Start key.  
The settings will be saved.  
The device switches off.

### Setting up a wireless group (Lrn)

To set up a wireless group, follow the instructions below.

1. Switch on the device.
2. Call up the menu.
3. Select the item “rF” from the menu.
4. Confirm your selection.

5. Select the “Lrn” (Learn) menu item from the “rF” submenu.
6. Confirm your selection.

The wireless group currently set (in this case: wireless group 0 “ld 0”) is displayed.

If wireless group “0” already exists, use the arrow keys to select a different ID (in this case: wireless group 1 “ld 1”).

7. Confirm your selection of wireless group.

The scale suggests a channel number for Channel 1 (in this case: “C1 0”).

You can adopt the proposed channel number or use the arrow keys to set a different channel number.

8. Confirm your selection for Channel 1.

The device suggests a channel number for Channel 2 (in this case: “C230”).

You can adopt the proposed channel number or use the arrow keys to set a different channel number.

#### NOTE

Double-digit channel numbers are displayed without a space. The display “C230” therefore means: Channel “2”, channel number “30”.

9. Confirm your selection for Channel 2.

The device suggests a channel number for Channel 3 (in this case: “C360”).

You can adopt the proposed channel number or use the arrow keys to set a different channel number.

10. Confirm your selection for Channel 3.

The message “StOP” appears in the display.

The device waits for signals from other wireless-capable devices in range.

#### NOTE

On some devices, a special switch-on procedure has to be followed if they are going to be integrated in a wireless group. Follow the instructions for use for the device in question.

11. Switch on the device, e.g. a wireless printer, you want to integrate in the wireless group.

When the wireless printer is detected, you will hear a beep.

#### NOTE

As soon as you have integrated a wireless printer in the wireless group, you then have to select a print option (Menu\rF\APrt) and set the time (Menu\rF\TIME).

12. Repeat step 11. for all the devices you want to integrate in this wireless group.
13. Press the Enter key to stop the search procedure.
14. Press one of the arrow keys to have the devices which have been detected displayed (in this case: "MO 3" for a wireless printer).

If you have integrated several devices in one wireless group, press the arrow keys several times to ensure that all devices have been detected by the scale.

15. Exit the menu using the Enter key or wait until you exit the menu automatically.
16. If you do not wish to make any more settings, press the Start key.  
The settings will be saved.  
The device switches off.

MO 3

### Activating automatic send (ASend)

You can configure the device so that the measurement results are automatically sent to all receivers which are ready to receive and are registered to the same wireless group (e. g.: wireless printer, PC with USB wireless adapter).

#### NOTE

If you are using a wireless printer, ensure that "Off" is not set as the print option (see "Selecting print option (APrt)" on page 105).

1. Switch on the device.
2. Select the "ASend" menu item from the "rF" submenu and confirm your selection.
3. Select the "On" setting and confirm your selection.  
You exit the menu automatically.

ASend

On

4. If you do not wish to make any more settings, press the Start key.  
The settings will be saved.  
The device switches off.

## Selecting print option (APrT)

You can configure the device so that the measurement results are automatically printed out on a wireless printer registered to the wireless group.

### NOTE

This function is only accessible if a seca wireless printer has been integrated in the wireless group via the “Learn” function.

1. Switch on the device.
2. Select the “APrT” menu item from the “rF” submenu and confirm your selection.
3. Select the setting relevant for your device combination:
  - HI: measurement results from length measuring devices
  - MA: measurement results from scales
  - HI\_MA: measurement results from length measuring devices and scales
  - Off: no automatic printing, printing only possible by pressing and holding the Enter key during the weighing procedure.
4. Confirm your selection.  
You exit the menu automatically.
5. If you do not wish to make any more settings, press the Start key.  
The settings will be saved.  
The device switches off.

## Setting the time (tIME)

You can configure the system so that the wireless printer automatically adds a date and time to your measurement results. To do this, you have to set date and time once on this device and transmit it to the internal clock of the wireless printer.

### NOTE

This function is only accessible if a seca wireless printer has been integrated in the wireless group via the “Learn” function.

1. Switch on the device.

2. Select the “tIME” menu item from the “rF” submenu.
3. Confirm your selection.  
The current setting for “Year” (**YEA**) is displayed.
4. Set the correct year.
5. Confirm your selection.
6. Repeat steps 3. and 4. accordingly for “month” (**Mon**), “day” (**dAy**), “hour” (**hour**) and “minute” (**Min**).
7. Confirm each of your selections.  
After confirming the setting for minute, you automatically exit the menu.  
The settings are automatically transmitted to the wireless printer.  
The wireless printer automatically adds a date and time to every printout.

#### NOTE

Follow the instructions for use for the wireless printer for further information about operating it.

8. If you do not wish to make any more settings, press the Start key.  
The settings will be saved.  
The device switches off.

## 7. HYGIENE TREATMENT

#### WARNING!

##### Electric shock

The device is not de-energized when the on/off key is pressed and the display goes out. Use of fluids on the device may cause electric shock.

- ▶ Ensure that the device is switched off before performing any hygiene treatment.
- ▶ Disconnect the power supply connector before performing any hygiene treatment.
- ▶ Before each hygiene treatment, take the rechargeable battery out of the device (if present and removable).
- ▶ Ensure that no fluids penetrate the device.

#### CAUTION!

##### Damage to device

Inappropriate detergents and disinfectants may damage the sensitive surfaces of the device.

- ▶ Use only disinfectants free of chlorine and alcohol which are explicitly suitable for acrylic sheet and other sensitive surfaces (active ingredient: quaternary ammonium compounds, for example).
- ▶ Do not use caustic or abrasive detergents.
- ▶ Do not use organic solvents (e.g. white spirit or petroleum spirit).

## 7.1 Cleaning

- ▶ Use a soft cloth dampened with mild soapsuds to clean the surfaces of the device.

## 7.2 Disinfecting

1. Check that your disinfectant is suitable for sensitive surfaces and acrylic sheet.
2. Follow the instructions on the disinfectant.
3. Disinfect the device:
  - ▶ Moisten a soft cloth with disinfectant and wipe down the device with it.
  - ▶ Comply with the intervals, see table.

| Interval                        | Component  |
|---------------------------------|--|
| <b>Before</b> every measurement | Weighing platform  |
| <b>After</b> every measurement  | Weighing platform  |
| If required                     | <ul style="list-style-type: none"> <li>• Display</li> <li>• Membrane keyboard</li> </ul> |

## 7.3 Sterilizing

This device may not be sterilised.

## 8. FUNCTION CHECK

- ▶ Perform a function check before each use.

A complete function check includes:

- a visual inspection for mechanical damage
- a visual inspection and function check of the display
- function check of all the controls shown in the section entitled “Overview” on page 71
- a function check of the optional accessories

If you find faults or deviations during the function check, as a first step try to remedy the fault with the aid of the section entitled “Troubleshooting” from page 108.

### **CAUTION!**

#### **Bodily injury**

If you find faults or deviations during the function check which you are unable to remedy with the aid of the section entitled “Troubleshooting” from page 108, you may not use the device.

- ▶ Have servicing and repairs carried out exclusively by seca service or by an authorized service partner.
- ▶ Please see the section entitled “Servicing/Verification” on page 110.

## 9. TROUBLESHOOTING

| <b>Problem</b>  | <b>Cause / Solution</b>  |
|---|--|
| <b>No weight displayed with a load on the scale.</b>  | The scale has no power supply.<br>- Check whether the scale is switched on   |
| <b>0.0 does not appear before weighing.</b>   | A load was already on the scale before it was switched on.<br>- Remove the load from the scale<br>- Switch the scale off and back on again |
| <b>One segment is illuminated either continuously or not at all.</b>  | The corresponding point has a fault.<br>- Inform the service department  |
| <b>The display  appears.</b> | Battery block power is running down.<br>- Recharge battery block as soon as possible   |
| <b>The display “bAtt” appears.</b>  | Battery block discharged.<br>- Charge battery block  |

| Problem  | Cause / Solution  |
|--|---|
| <b>The display “STOP” appears.</b>   | Maximum load has been exceeded.<br>- Remove the load from the scale   |
| <b>The display “tEMP” appears.</b>   | The ambient temperature of the scale is too high or too low.<br>- Set up the scale in an ambient temperature between +10 °C and +40 °C<br>- Wait for around 15 minutes until the device has adapted to the ambient temperature  |
| <b>After switching on for the first time, measurement results are transmitted and two beeps are heard.</b> | <ul style="list-style-type: none"> <li>• The device was unable to send measurement results to the wireless receiver (seca wireless printer or PC with seca USB wireless adapter).<br/>- Ensure that the scale is integrated in the wireless network<br/>- Ensure that the receiver is switched on</li> <li>• Reception is being interfered with by HF equipment (e.g. cell phones) in the vicinity.<br/>- Keep HF equipment at least 1 meter away from transmitters and receivers in the seca wireless network. The actual transmission output of HF equipment may require minimum distances of more than 1 meter. For details, go to <a href="http://www.seca.com">www.seca.com</a>.</li> </ul> <p><b>NOTE</b><br/>If this fault is not eliminated, there will be no fresh acoustic warning when further attempts are made to send data.</p> |
| <b>Only the “SYS” item is visible in the “rF” menu.</b>  | The wireless module is deactivated.<br>- Activate wireless module (see “Activating the wireless module (SYS)” on page 102)  |
| <b>Only the “SYS” and “Lrn” items are visible in the “rF” menu.</b>  | The wireless module is activated and no wireless group is set up.<br>- Set up wireless group (see “Setting up a wireless group (Lrn)” on page 102)  |
| <b>The items “APrt” and “tIME” are not visible in the “rF” menu.</b>                                       | No wireless printer registered to the wireless group.<br>- Register the wireless printer to the wireless group via the “Lrn” menu item (see “Setting up a wireless group (Lrn)” on page 102)  |
| <b>After calling up the menu, the “rF” item is not displayed.</b>  | The scale’s wireless module is defective.<br>- Inform seca service  |

| Problem   | Cause / Solution   |
|---|--|
| <b>The display “Er:[No.]:11” appears.</b>   | The scale has too high a load or too high a load in one corner.<br><ul style="list-style-type: none"> <li>- Take the load off the scale or distribute the weight more evenly</li> <li>- Restart the scale</li> </ul> |
| <b>The display “Er:[No.]:12” appears.</b>   | The scale has been switched on with too high a load.<br><ul style="list-style-type: none"> <li>- Remove the load from the scale</li> <li>- Restart the scale</li> </ul>  |
| <b>The display “Er:[No.]:16” appears.</b>   | The scale was caused to oscillate, the zero point could not be determined.<br><ul style="list-style-type: none"> <li>- Restart the scale</li> </ul>  |
| <b>The Enter key (send/print) is pressed and the display “Er:[No.]:71” appears.</b> | No data transmission possible, wireless module deactivated.<br><ul style="list-style-type: none"> <li>- Activate wireless module (see “Activating the wireless module (SYS)” on page 102)</li> </ul>                 |
| <b>The Enter key (send/print) is pressed and the display “Er:[No.]:72” appears.</b> | No data transmission possible, no wireless group set up.<br><ul style="list-style-type: none"> <li>- Set up wireless group (see “Setting up a wireless group (Lrn)” on page 102)</li> </ul>                          |

## 10.SERVICING/VERIFICATION

### 10.1 Information on servicing and verification

We recommend having your device serviced prior to verification.

**NOTICE!**

**Incorrect measurements as a result of poor servicing**

- ▶ Have servicing and repairs carried out exclusively by seca Service or by an authorized service partner.
- ▶ You can find service partners in your area at [www.seca.com](http://www.seca.com) or by sending an e-mail to [service@seca.com](mailto:service@seca.com).

Have an authorized technician perform verification according to national legal regulations.

Verification is necessary whenever one or more verification seals are damaged or the contents of the verification counter no longer match the number on the applicable verification counter sticker. If verification seals are damaged, contact seca Service directly.

## 10.2 Checking the content of the verification counter


This seca scale is verified. Verifications may only be performed by authorized agencies. To guarantee this, the scale is equipped with a verification counter to record any change in verification-related data.

If you want to check whether the scale has been properly verified, proceed as follows:

1. If necessary, turn off the scale.
2. Press and hold any key and start the scale.

The current content of the verification counter flashes in the display for a few seconds.

3. Compare the content of the verification counter output with the number given on the verification counter sticker.

Both numbers have to match for the verification to be valid. If the sticker and the verification counter do not match, the scale must be verified. Please contact your service partner or seca Service. Once the scale has been verified, a new, updated verification counter sticker is used to identify the status of the verification counter. The person authorized to perform the verification secures this sticker with an additional seal. The verification counter sticker can be ordered from seca Service.

# 11. TECHNICAL DATA

## 11.1 General technical data

| General technical data seca 645   |   |
|---|---|
| Dimensions of scale <ul style="list-style-type: none"> <li>• Depth</li> <li>• Width</li> <li>• Height</li> </ul>  | 750 mm<br>750 mm<br>1140 mm   |
| Dimensions of weighing platform <ul style="list-style-type: none"> <li>• Depth</li> <li>• Width</li> <li>• Height</li> </ul>  | 600 mm<br>600 mm<br>45 mm   |
| Net weight  | approx. 24 kg   |
| Temperature range <ul style="list-style-type: none"> <li>• Operation</li> <li>• Storage</li> <li>• Transport</li> </ul>   | +10 °C to +40 °C / +50 °F to 104 °F<br>-10 °C to +65 °C / +14 °F to 149 °F<br>-10 °C to +65 °C / +14 °F to 149 °F |
| Air pressure <ul style="list-style-type: none"> <li>• Operation</li> <li>• Storage</li> <li>• Transport</li> </ul>  | 700-1060 hPa<br>700-1060 hPa<br>700-1060 hPa  |
| Humidity <ul style="list-style-type: none"> <li>• Operation</li> <li>• Storage</li> <li>• Transport</li> </ul>  | 30 % - 80 % no condensation<br>0 % - 95 % no condensation<br>0 % - 95 % no condensation                           |
| Height of digits  | 25 mm   |
| Power supply <ul style="list-style-type: none"> <li>• Battery block               <ul style="list-style-type: none"> <li>- Supply voltage</li> <li>- Maximum power consumption</li> </ul> </li> <li>• Power pack (model-dependent)               <ul style="list-style-type: none"> <li>- Supply voltage</li> <li>- Power consumption with wireless module deactivated</li> <li>- Power consumption with wireless module activated</li> <li>- Insulated device, protection class II (EN 60601-1)</li> </ul> </li> </ul> | 7.2 V<br>Typically 120 mA<br>12 V<br>Approx. 25 mA<br>Approx. 42 mA<br><input type="checkbox"/>                   |
| Maximum runtime in battery mode <ul style="list-style-type: none"> <li>• with wireless module deactivated</li> <li>• with wireless module activated</li> </ul>  | approx. 2,400 weighing operations<br>Power pack recommended   |

| General technical data seca 645   |   |
|---|---|
| Medical device in accordance with Directive 93/42/EEC   | Class I with measuring function   |
| EN 60601-1: electrical medical device, type B   |  |
| Wireless transmission <ul style="list-style-type: none"> <li>• Frequency band</li> <li>• Transmission power</li> <li>• Standards applied</li> </ul> | 2.433 GHz - 2.480 GHz<br>< 10 mW<br>EN 300 328<br>EN 301 489-1<br>EN 301 489-17   |

## 11.2 Technical data, weighing

| Technical data, weighing seca 645  |                                     |
|--|-------------------------------------|
| Calibration in accordance with Directive 2014/31/EU  | Class III                           |
| Maximum load <ul style="list-style-type: none"> <li>• Partial weighing range 1</li> <li>• Partial weighing range 2</li> </ul>  | 200 kg<br>300 kg                    |
| Minimum load <ul style="list-style-type: none"> <li>• Partial weighing range 1</li> <li>• Partial weighing range 2</li> </ul>  | 2.0 kg<br>4.0 kg                    |
| Increments <ul style="list-style-type: none"> <li>• Partial weighing range 1</li> <li>• Partial weighing range 2</li> </ul>  | 100 g<br>200 g                      |
| Tare range   | 300 kg                              |
| Accuracy on initial calibration <ul style="list-style-type: none"> <li>• Partial weighing range 1, to 50 kg</li> <li>• Partial weighing range 1, 50 kg to 200 kg</li> <li>• Partial weighing range 2, to 100 kg</li> <li>• Partial weighing range 2, 100 kg to 300 kg</li> </ul> | ±50 g<br>±100 g<br>±100 g<br>±200 g |

## 12. OPTIONAL ACCESSORIES

| seca 360° wireless devices  | Article number   |
|---|--|
| Length measuring rods <ul style="list-style-type: none"> <li>• <b>seca 274</b></li> <li>• <b>seca 264</b></li> </ul>  | Country-specific versions<br>Country-specific versions |
| Wireless printers <ul style="list-style-type: none"> <li>• <b>seca 360° wireless printer 465</b></li> <li>• <b>seca 360° wireless printer advanced 466</b></li> </ul> | Country-specific versions<br>Country-specific versions |

| <b>seca 360° wireless devices</b>          | <b>Article number</b>                 |
|--|---------------------------------------|
| PC software<br>• <b>seca analytics 115</b> | Application-specific license packages |
| <b>seca 360° wireless USB adapter 456</b>  | 456-00-00-009                         |

## 13.SPARE PARTS

| <b>Spare parts</b>  | <b>Article number</b> |
|---|-----------------------|
| Power pack with Euro-connector: 230 V~ /50 Hz/12 V= / 130 mA                  | 68-32-10-252          |
| Switchmode power pack with adapters: 100 - 240 V~ / 50 -60 Hz / 12 V= / 0.5 A | 68-32-10-265          |

## 14.DISPOSAL

### 14.1 Device

Do not dispose of the device with household waste. The device must be disposed of properly as electronic waste. Comply with the national provisions applicable in your country. For further information contact our service department at:

**service@seca.com**

### 14.2 Batteries and rechargeable batteries

Spent (rechargeable) batteries should not be discarded with household waste, regardless of whether they contain harmful substances or not. As a consumer you are obliged by law to dispose of (rechargeable) batteries via the collection points set up by the municipal authorities or the retail sector. Only discard (rechargeable) batteries when fully discharged.

## 15. WARRANTY

We offer a two-year warranty from the date of delivery for defects attributable to faulty material or poor workmanship. This excludes all moveable parts such as (rechargeable) batteries, cables, power supply units, etc. Defects which are covered by the warranty shall be rectified free of charge for customers on production of the sales receipt. No further claims can be accepted. The costs of shipment in both directions shall be borne by the customer where the device is not located at the customer's premises. In the event of any damage during shipment warranty claims can only be asserted where the complete original packaging was used for shipment and the scales were secured inside in the same manner as in the original packaging. You should therefore keep all packaging.

The warranty shall become null and void where the device is opened by persons not expressly authorised to do so by seca.

In the event of a warranty issue, please contact your local seca office or the dealer from whom you ordered the product.

## 16. DECLARATION OF CONFORMITY

seca gmbh & co. kg hereby declares that the product meets the terms of the applicable European directives. The unabridged declaration of conformity can be found at: [www.seca.com](http://www.seca.com).

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A rectangular digital display showing the text "AHOLD" in a seven-segment font.

A rectangular digital display showing the text "On" in a seven-segment font.

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A rectangular digital display showing the text "bEEP" in a seven-segment font.

A rectangular digital display showing the text "PrESS" in a seven-segment font.

A rectangular digital display showing the text "On" in a seven-segment font.

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# DECLARATION OF CONFORMITY FOR USA AND CANADA

seca  
seca 645

FCC ID: X6T172A01

IC: 8898A-172A01

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions. (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

## NOTE

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

## NOTE

Changes or modifications made to this equipment not expressly approved by seca may void the FCC authorization to operate this equipment.

## NOTE

Radiofrequency radiation exposure Information: This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance of 1 m between the radiator and your body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

# Medical Measuring Systems and Scales since 1840

seca gmbh & co. kg  
Hammer Steindamm 3–25  
22089 Hamburg · Germany  
Telephone +49 40 20 00 00 0  
Fax +49 40 20 00 00 50  
info@seca.com

seca operates worldwide with headquarters  
in Germany and branches in:

**seca** france

**seca** united kingdom

**seca** north america

**seca** schweiz

**seca** zhong guo

**seca** nihon

**seca** mexico

**seca** austria

**seca** polska

**seca** middle east

**seca** brasil

**seca** suomi

**seca** américa latina

**seca** asia pacific

and with exclusive partners in  
more than 110 countries.

All contact data under [www.seca.com](http://www.seca.com)

**seca**<sup>®</sup>  
Precision for health